





Ph: 888-777-8835 • 209-839-2990 • Fax: 209-839-2999
sales@m-depot.com • www.microscope-depot.com
392 W. Larch Rd., Bldg. 28, Tracy, CA. 95304 USA



S-14700 - 1" Total Scale • **S-14702** - same in Reflected Light version

0 thru .9 - divided into 0.100" divisions / .9 thru 1.0 divided into 0.010" divisions / last two 0.010" divisions divided into 0.001" divisions



S-14710 - 25mm Stage Micrometer / 1.0mm - .1mm - .02mm Div. • **S-14712** - same in Reflected Light version

0 thru 20 divided into 1.0mm divisions / 20 thru 25 divided into 0.1mm divisions / 23 thru 25 further divided into 0.02mm divisions

S-14720 - 1" & 25mm Stage Micrometer / both above scales on same slide • **S-14722** - same in Reflected Light version

SPECS: Glass: 1.0" x 3.0" x 0.110" / Lines: Chrome deposition / Lines on RF version: Clear on chrome background
Flatness: 2 waves or better / Tolerances: ± 0.005 mm unless otherwise noted / Surface Quality: 20/10 / Parallelism: 2 min.
Calibration traceable to NIST available at additional cost