



NATIONAL MODEL 186, 187, 188 & 189 SHOP MICROSCOPE SERIES

See illustration on page 2.

Compact design and precision optics combine in this indispensable shop tool. Small size and rugged carrying case make it easy to take anywhere on the job. Four models provide a choice of four different magnifications. Each has the following features:

OPTICAL SYSTEM

- WF10x eyepiece with direct reading reticle (no calibration required). Eyepiece turns to permit sharp focus of reticle scale.
- Each model includes both inches and metric scales.
- See chart below for magnification and other optical specifications.

ILLUMINATION

- Penlight illuminator in swivel mount. Requires 2 AA batteries, included.

FOCUSING

- Rack & pinion focusing with dual focusing knobs.

FRAME

- Cast aluminum, gray enamel finish.
- Height 230mm (9 in.). Net wt. .64 kgs (1.4 lbs).

PACKAGING

- Rugged, molded black plastic case with latch and carrying handles. Interior is fitted to securely accommodate microscope in place.
- Includes instruction brochure, dustcover and 2 AA batteries.
- Individual: Ship wt. 2.3 kgs (5 lbs),
- Master carton pack 10 each: Ship wt. 16.8 kgs (37 lbs),

WARRANTY

- 5-year limited warranty, excluding penlight and batteries.

MODEL NO.	MAGNIFICATION	METRIC SCALE	INCH SCALE	FIELD DIAMETER	WORKING DISTANCE
186	20x	0.06mm x 100	0.003" x 100	9.0mm	23mm
187	40x	0.03mm x 100	0.0015" x 100	4.5mm	23mm
188	60x	0.02mm x 100	0.001" x 100	3.0mm	23mm
189	100x	0.01mm x 100	0.0005" x 100	1.8mm	12.0mm

PARTS & ACCESSORIES:

610-186	WF10x Eyepiece with reticle, for Model 186
610-187	WF10x Eyepiece with reticle, for Model 187
610-188	WF10x Eyepiece with reticle, for Model 188
610-189	WF10x Eyepiece with reticle, for Model 189
610-190	WF10x Eyepiece without reticle
702-186	2x objective for Model 186
704-187	4x objective for Model 187
706-188	6x objective for Model 188
710-189	10x objective for Model 189
904	Penlight illuminator (batteries not included in replacement penlight)



NATIONAL MODEL 187

Distributed by:

the
microscope  *depot*

Sales: 888-777-8835 • Email: sales@m-depot.com

Main: 209-839-2990 • Fax: 209-839-2999

<http://www.microscope-depot.com>