



NATIONAL MODEL D-CAM-1 ANALOG CAMERA

This exciting and very affordable teaching tool is great for demonstrating microscopy to an entire class. The D-CAM-1 mounts on any microscope eyepiece and supplies a superior analog image to your video monitor through an RCA cable. This complete package ready for use comes with CS camera adapter, RCA cable; 10mm CCD lens; 4 eyepiece adapters that fit over microscope eyepieces with O.D. of 28mm, 30mm, 34mm and 35mm; 110-220v power converter and cord, Macro viewing tube, 1 year limited warranty.

Warranty: 1 Year Limited Warranty

Includes:

- Analog Camera, 1/3 in. Interline CCD
 - Analog signal 330 TV lines of resolution
 - Minimum illumination: 2 Lux
- CS camera mount.
- 12mm CCD lens.
- 4 ea. Eyepieces adaptors to fit eyepiece O.D.'s of 28mm, 30mm, 34mm, 35mm.
- RCA cable.
- 110v-240v power converter and power cord.
- Macro Tube

Specifications:

- Camera (NTSC system) has 1/3" CCD chip (330 lines TV resolution)
- Power supply: 12V power cord.
- Magnification: 50:1 and greater
- Power consumption: 2.28 W or less



NATIONAL MODEL #D-CAM-1



MODEL #D-CAM-1 – In Use (Microscope & TV not included)

Distributed by:

the
microscope  *depot*

Sales: 888-777-8835 • Email: sales@m-depot.com

Main: 209-839-2990 • Fax: 209-839-2999

<http://www.microscope-depot.com>