



MOTICAM TABLET

the next generation of digital imaging

Introducing a new line of MotiCam tablets designed to fit increasing demands of technology and digital imaging. With preloaded Moti apps, it is easy for teachers, students, and researchers to connect, view, and share live images. Using a C-mount adapter, these interactive tablets can be used with Swift Optical's **M10T**, **M15T**, **M3-F**, and **M29TZ** microscopes; and National Optical's **163** and **420T** microscopes.

SPECIFICATIONS:

Display	7" LCD Screen with 1200x800 *16:9 pixels resolution 10" LCD Screen with
Connection	WiFi; mini HDMI; micro SD card
Sensor	5.0 MP
Captured Image Resolution	5.0 MP
Captured Video Resolution	720P
HDMI	1080
Operating System	Android
WiFi	WiFi output (802.11 ab/g/n supporting both hot-spot mode or using existing WiFi network)
App	Preloaded optimized MotiConnect Supplied with Moti Images Plus Software for PC/Mac Built-in Moti Hub allows wireless sharing of MotiCam S2 and T2 real time images to multiple devices

